

Drug of Abuse Tests

Colloidal Gold Rapid Diagnostic Tests Qualitative Tests



FEATURES

- ✓ Detect upto 25 drugs and metabolites in urine.
- ✓ One step operation, results in 5 minutes.
- ✓ Certificated with CE, FDA 510k, CLIA waived.
- ✓ Wide arrange of applications.
- ✓ Room temperature (2~30°C) storage up to 24 months.
- ✓ 6 adulterants available: PH, SG, OX, CR, NI and Glutaraldehyde.



Test Item	Sensitivity (ng/ml)	Format	Specimen
ALC Alcohol	>0.02% BAC	Strip	Saliva
AMP Amphetamine	500/ 1000	Strip Single Dip Card Multi Panel Urine Cup	Urine
BAR Barbiturates	200/ 300		
BUP Buprenorphine	10		
BZO Benzodiazepines	300		
COC Cocaine	150/ 200/ 300		
MDMA Methylendioxyamphetamine Ecstasy	500		
MET Methamphetamine	500/ 1000		
MOP Morphine	100/ 200/ 300		
MTD Methadone	300		
OXY Oxycodone	100		
PCP Phencyclidine	25		
PPX Propoxyphene	300		
TCATri-cyclic Antidepressants	1,000		
THC Marijuana	50		
OPIOpiate 2000	2,000		
COTCotinine	100/ 200		
EtG Ethyl Glucuronide	500		
FYL Fentanyl	200		
K2 Synthetic Cannabinoid	50		
KET Ketamine	1,000		
6-MAM 6-Aoetylmorphine	10		
MQL Methaqualone	300		
TRATramadol	200		

More new test items coming based on market needs.

