

Tumor Marker Tests

Colloidal Gold Rapid Diagnostic Tests

Qualitative Tests



FEATURES

- ✓ High sensitivity and specificity.
- ✓ One step operation, results in 5 minutes.
- ✓ Wide arrange of applications.
- ✓ Room temperature (2~30°C) storage up to 36 months.

Test Item	Sensitivity (ng/ mL)	Format	Format
FOB	100/200	Strip Cassette	Whole blood/Serum/Plasma
PSA	4	Strip Cassette	Whole blood/Serum/Plasma
AFP	25	Strip Cassette	Whole blood/Serum/Plasma
CEA	5	Strip Cassette	Whole blood/Serum/Plasma