



FEATURES

- ✓ High sensitivity and specificity.
- ✓ One step operation, results in 5 minutes
- ✓ Wide arrange of applications.
- ✓ Room temperature [2~30°C]storage upto 24 months

Dog

CDV Ag *Canine Distemper Virus Antigen*

CDV Ab *Canine Distemper Virus Antibody*

CPV Ag *Canine Parvo Virus Antigen*

CPV Ab *Canine Parvovirus Antibody*

CIV Ag *Canine Influenza Virus Antigen*

CCV Ag *Canine Coronavirus Antigen*

CAV Ag *Canine Infectious Hepatitis Antigen*

CDV/CAV Ag *Canine Distemper/Adeno Antigen*

CPV/CDV Ab *CPV/CDV Antibody*

CPV/CDV Ag *CPV/CDV Antigen*

CPV/CCV Ag *Canine Parvo/Corona Antigen*

C.Prog *Canine Progesterone*

C.CRP *Canine CRP*

Giardia Ag *Canine Giardia Antigen*

CPV/CCV/Giardia Ag *CPV/CCV/Giardia Antigen*

CDV/CAV/CIV Ag *Canine Distemper/Adeno/Influenza Antigen*

CHW Ag *Canine Heartworm Antigen*

E.canis Ab *Ehrlichia Canis Antibody*

BG Ab *Babesia Gibsoni Antibody*

Lyme Ab *Lyme Disease Antibody*

RLN *Canine Pregnancy Relaxin*

Rabies Ab *Rabies Antibody*

Rabies Ag *Rabies Antigen*

C.Bru Ab *Canine Brucella Antibody*

LSH Ab *Leishmania Canis Antibody*

Anaplasma Ab *Anaplasma Antibody*

EHR/ANA *Ehrlichia/Anaplasma*

EHR/BAB/ANA *Ehrlichia/Babesia/Anaplasma*

EHR/LYM/ANA *Ehrlichia/Lyme/Anaplasma*

EHR/BAB/ANA/CHW *Ehrlichia/Babesia/Anaplasma/Heartorm*

EHR/LYM/ANA/CHW *Ehrlichia/Lyme/Anaplasma/Heartorm*

EHR/LSH/ANA/CHW *Ehrlichia/Leishmania/Anaplasma/Heartorm*